

PRODIGY[®]

Count-a-dose[®]

The Simple Way to Fill Insulin Syringes

Features

The Prodigy[®] Count-a-dose™ is a simple way to fill insulin syringes safely and accurately with one or two bottles of insulin (any brand). The bottle holder is marked to identify bottles for easy mixing—one raised dot for the first bottle and two raised dots for the second bottle. With the click dial starting all the way down toward the raised minus (-) sign, the user simply moves the click dial the number of times, or clicks, per unit of insulin needed. Each click of the dial toward the plus sign (+) will draw one (1) unit of insulin into the syringe.

- Low vision and blind friendly
- Measures one (1) unit of insulin at a time, allowing better control of insulin therapy
- Makes a distinctive “click” with each increment
- Holds one or two bottles of insulin—any brand for easy mixing
- Lightweight; easy to carry anywhere



Fill Syringes with Independence

The Prodigy[®] Count-a-dose™ is a simple system to learn, is well-adapted to the use of non-sighted alternative techniques and allows a blind diabetic to measure from one to fifty units of insulin in single-unit increments. The Prodigy[®] Count-a-dose™ uses a Prodigy[®] 50-unit (1/2 cc) syringe or a BD (Becton Dickinson) 50-unit (1/2 cc) insulin syringe.

www.prodigymeter.com/countadose

PRODIGY[®]

Count-a-dose[®]

The Simple Way to Fill Insulin Syringes

+ Notch

Tactile raised plus (+) sign indicates the direction to move click dial to increase amount of insulin in single unit increments.

Click Dial

Increases or decreases the amount of insulin in single-unit increments.

- Notch

Tactile raised minus (-) sign indicates the direction to move click dial to decrease amount of insulin in single unit increments.

T-Shaped Bar

Holds two bottles of insulin and syringe needle.

Insulin Identification Dots


Raised dots located on both sides of the T-Bar used to distinguish between the two bottles of insulin.

Syringe Clip

Secures the barrel of the syringe.

Flange Notch

Secures the flange of the syringe.



To learn more about Prodigy[®] Count-a-dose[®], call us toll-free at **1.866.908.9201** or visit us online at **www.prodigymeter.com**.

PRODIGY[®]

9300 Harris Corners Parkway, Suite 450
Charlotte, North Carolina 28269