

Prodigy Voice™ Blood Glucose Monitoring System

Features

The Prodigy Voice™ totally audible blood glucose meter is excellent for the blind and low vision diabetic patient. The Prodigy Voice™ requires no coding, which means fewer steps and avoids the danger of false results caused by coding errors. Simply insert the test strip into the strip port—located at the top of the meter—and the Prodigy Voice™ automatically turns on. It talks the user through the testing steps, making it easy to use while providing safe, accurate test results every time.

- Easy to use—no coding required
- Totally audible
- Approved for alternate site testing
- Hear and see accurate results in just 7 seconds
- Standard USB port to easily download test results


Award Winning Meter

The Prodigy Voice™ meter is made with quality components, comes with a lifetime warranty and is clinically proven to have higher accuracy than standard FDA requirements.

The Prodigy Voice™ meter is the only meter to receive Access Awards from both the National Federation of the Blind (NFB) and American Foundation for the Blind (AFB). The term Access Plus means that the product is not just accessible to blind people, but also useful and convenient. “The A+ is for a product that affords the blind the same convenience and features available to everybody else.”¹

The Prodigy Voice™ was carefully designed with help from the NFB as well as Certified Diabetes Educators. It talks the user through setup steps, testing results and memory recall. Its tactile features allow the user to easily feel and locate the different buttons...allowing complete independence for the blind or low vision patient.

www.prodigymeter.com/voice



Prodigy Voice™

Totally Audible...Complete Independence!



AFB American Foundation™
for the Blind

Test Strip Port

The indentation makes it easy for fingers to locate where to insert a test strip. Then the meter automatically turns on.

Time and Date

Set Button

Allows the user to setup the functions of the meter using audible prompts.

Repeat Button

Allows the user to replay the last message or test result.

Ear Phone Jack

Allows the user to test with audible privacy.

Speaker

A clear, loud voice is heard.

Mini USB Port

Uses any standard USB Cable. Download test results to a computer with Prodigy® meter's FREE software.

To learn more about the Prodigy Voice™, call us toll-free at **1.800.366.5901**
or visit us online at **www.prodigymeter.com**

PRODIGY®

9300 Harris Corners Parkway, Suite 450
Charlotte, North Carolina 28269